

Enjoy life's true sound.

Equipped with the next generation of BestSound™ Technology, Motion provides binaural directional performance that is so advanced, wearers can understand speech in demanding environments even better than people with normal hearing.* This world-first technology delivers rich, clear and natural sound quality.

Hearing with two ears, or binaural hearing, helps us to localize sound sources, maintain spatial orientation, and understand speech better. Inspired by these advantages, Motion's cutting-edge technology delivers an exceptional binaural listening experience.

Siemens Motion models do all this with the industry's lowest battery consumption.** This provides the best listening experience, all the time, without concerns about battery life.

* Two independent clinical studies have shown that Siemens Pure® and Carat™ hearing instruments with binaural processing provide better than normal hearing capabilities in certain demanding listening environments. The new wireless, two microphone Motion has the same directional performance as the Pure and Carat binax.

** For instruments with binaural processing.



iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

Android and Google Play are trademarks of Google Inc. Apple App Store is a trademark of Apple Inc.

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos GmbH is under license. Other trademarks and trade names are those of their respective owners.

The information in this brochure is subject to change without prior notice. © 2015 Sivantos GmbH. All rights reserved.

Hearing instruments help many people hear better, but no hearing instrument can solve every hearing problem nor restore normal hearing.

Legal Manufacturer

Sivantos GmbH
Henri-Dunant-Strasse 100
DE-91058 Erlangen
Germany

Sivantos GmbH is a Trademark Licensee of Siemens AG.

SIEMENS

Motion



usa.siemens.com/motion

Enjoy the sound of life.

The new Motion® hearing aids.
Versatile, comfortable, and exceptional binaural hearing.



Sleek, easy to handle and exceptional, binaural hearing:

The new Motion.


Equipped with innovative sleek housings, all new Motion® models discreetly hide behind your ear and offer maximum wearing and handling comfort. The smart design also enhances directional speech understanding. This way, Motion lets you understand your conversations – even in demanding listening situations such as busy restaurants.



Local Contact

Sivantos, Inc.
10 Constitution Avenue
Piscataway, NJ 08854
www.sivantos.com

Read the QR code with the QR code reader on your smart phone



usa.siemens.com/hearing

SH201503 2/15 1.0 22B771 2K SHI/15702-15



Life sounds brilliant.

Reliable, hassle-free and versatile:
Motion BTEs.



As behind-the-ear hearing aids, Motion models are very robust. IP67-rated, they are resistant to moisture, sweat, dust, and dirt – for greater reliability.



Exceptional binaural hearing –
in any situation:
The new Motion.

No matter what your day has in store, you don't want to think about your hearing aids. All Motion models offer a perfect solution to almost all hearing needs and preferences. The rechargeable Motion SX and Motion SA with optional audio shoe for FM systems, come with a new **sleek housing** that completely hides behind your ear.

And the rechargeable Motion PX offers **high amplification**. This makes it a perfect solution even for wearers with severe hearing loss. Plus, no more changing tiny batteries – Motion PX and Motion SX can be recharged, refreshed, and dehumidified overnight with the optional Siemens eCharger™.

All Motion models include a telecoil for better hearing in facilities with an induction loop, such as places of worship, theaters, or public buildings.

Hearing that is more precise and clear.

Immerse yourself in conversations without a second thought because Motion works automatically – even in environments with a lot of noise around you. Motion focuses precisely on speech coming from in front of you, so you can enjoy every word effortlessly. Or, if you're out for a walk and somebody is

talking to you from the side, Motion automatically adjusts so you can take full pleasure in your conversation. It even helps you enjoy every minute outdoors by preventing annoying wind noise.

Getting the better of tinnitus.

Tinnitus can be distracting and irritating, especially in quiet environments. With Motion's optional **tinnitus therapy signals**, relief can be provided by a personalized therapy signal that can reduce the effects of tinnitus, thereby bringing back comfort and serenity.

The ultimate in seamless
connections and discreet control.

touchControl App.

Now, all Motion binax hearing aids can be controlled with a smartphone. The touchControl™ App turns an Android™ or Apple® iOS® smartphone into a hearing aid remote control. Adjusting volume, treble, and selecting programs is now possible with a FREE app – no other intermediary device is needed.



touchControl App


easyTek remote audio streamer.

With a one-button press, the easyTek™ remote audio streamer transmits audio from your TV, phone and *Bluetooth*®-enabled devices right into your hearing aids. Worn around the neck or hidden under clothing, the lightweight easyTek transforms Motion wireless hearing aids into a high-quality stereo headset, pairing with your favorite devices right out of the box.



easyTek App.

Combined with the easyTek, the easyTek™ App turns a smartphone into a hearing aid control center enabling switching between audio sources or 'zooming-in' on a conversation at a noisy restaurant – it's as easy and discreet as sending a text.



easyTek

easyTek App

The easyTek is compatible with Android, Apple iOS and Windows devices.

Use your smartphone to control
how and what you want to hear.

Even though your hearing aids automatically adjust to your listening environments, with easyTek and the easyTek App, you now have the ability to manually adjust them to suit your unique listening situation with your smartphone app.

Spatial Configurator Span.

Using the slider at the bottom of the easyTek App, you can slide right to focus the sound from the front while suppressing the sound from other directions. Or, you can slide to the left to hear the entire surrounding environment.



Spatial Configurator Direction.

You can adjust the hearing aids to focus on the sound or person you want to hear, including front, back, or either side of you. Just press the direction from which you want to hear to activate Spatial Configurator Direction.

You can also use the Spatial Configurator via a rocker switch on your hearing aids. Your hearing care professional can show you how.



Additional accessories.

The easyPocket™ remote control and the VoiceLink™ companion microphone further enhance Motion's easy handling. There is an optional audio shoe for connecting to FM systems.



Audio shoe

VoiceLink

Motion at a glance:



Motion's range of colors are designed to match a variety of hair and skin tones.

| | | | |
|---------------|----------------|--------------|----------------|
| Beige | Granite | Grey | Black |
| Silver | Sandy Brown | Dark Granite | Brown |
| Golden Blonde | Dark Champagne | Candy Pink* | Galactic Blue* |
| Elegance* | Pearl White* | Spirit* | |

*Additional colors for some models.

